ENGLISH
LIGHT-CURE FLUORIDE RELEASING DENTAL FLOWABLE RESTORATIVE MATERIAL

Description
BEAUTIFIL Flow is a radiopaque light-cure fluoride releasing flowable restorative material for esthetic restoration containing S-PRG filler created by using the Pre-Reacted Glass-ionomer (PRG) technology. BEAUTIFIL Flow is one of products of "GIOMER", a product category that hybridized the advantages of glass-ionomer and composite resin. BEAUTIFIL Flow consists of two types, F02 with low flow and F10 with high flow, which can be utilized for different clinical case. BEAUTIFIL Flow (F02 & F10) has excellent handling based on flow properties and a novel syringe and excellent esthetics by controlling light diffusion.

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Indications

- Minimally invasive cavity restorations
- Anterior restorations such as class III and class V involving cervical caries, root surface caries, wedge-shaped defects
- Small posterior restorations such as class I and class II without occlusal stress bearing
- Deciduous restorations
- Developmental teeth defects like enamel hypoplasia restorations
- Preventive restoration (pit & fissure sealant)
- Base/liner under restoration
- Undercure blockout
- Additional or marginal filling and recontouring after composite resin restoration
- Color adjustment on the surface of composite resin restoration and discolored teeth
- Direct composite resin laminate veneer
- Cementation of laminate veneer

Shades
A1, A2, A3, A3.5, A4, A3O
G (Gum shade) for F02 for cervical area
A3T (A3 Translucency) for F10

Clinical depth of cure			
Shade	A2, A3T	A1, A3, G	A3.5, A4, A3O
Clinical depth of cure (in mm)	3,0	2,5	2,0
For the above value is measured with using a tungsten-halogen lamp polymerization activator for 30 sec.			

- For Direct filling**
- After cleaning and removing the stains from the teeth, take shade while the tooth is still moist.
 - Make cavity preparations in conventional method. If a pulp exposure occurred or if a cavity is deep, cover the areas with a calcium hydroxide material. (Avoid using eugenol-containing liner.)
 - To promote adhesion the use of SHO-FU's FL-Bond is preferably recommended for double action of fluoride release.
 - In any case follow manufacturer's instructions.
 - Pull the shield cap off from a syringe and attach a needle tip by turning it clockwise until it is securely attached. Dispense necessary amount of BEAUTIFIL Flow directly into the cavity. BEAUTIFIL Flow may also be extruded onto a dispensing pad to be applied into cavity with a suitable instrument. Incremental filling of 2mm is recommended. Remove the needle tip from a syringe by turning it counter-clockwise and discard. Then, attach a new needle tip and reattach the opposite side of shield cap to the tip of a needle tip.
 - After applying, light cure each layer entirely in fanning motion for 30 seconds using a curing unit (halogen lamp) or for 10 seconds using a high intensity visible light with minimum intensity of 1000 mW/cm2. Light curing time may differ depending on shade and depth of the cavity (refer to the above chart for the depth of cure and shade). When the cavity surface is large, light cure in segment.
 - Remove excess material with suitable finishers such as "Dura White Stones", "T&F Hybrid Points", or "Super-Snap Strips" for interproximal removal and after adjusting occlusion, use polishers such as "CompoMaster", "Super-Snaps" or "OneGloss" for final polishing.

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- For Base/lining**
- After cleaning and removing the stains from the teeth, take shade while the tooth is still moist.
 - Make cavity preparations in conventional method. If a pulp exposure occurred or if a cavity is deep, cover the areas with a calcium hydroxide material. (Avoid using eugenol-containing liner.)
 - To promote adhesion the use of SHO-FU's FL-Bond is preferably recommended for double action of fluoride release.
 - In any case follow manufacturer's instructions.
 - Pull the shield cap off from a syringe and attach a needle tip by turning it dckwise until it is securely attached. Dispense necessary amount of BEAUTIFIL Flow directly into the cavity. BEAUTIFIL Flow may also be extruded onto a dispensing pad to be applied into cavity with a suitable instrument. Incremental filling of 2mm is recommended. Remove the needle tip from a syringe by turning it counter-clockwise and discard. Then, attach a new needle tip and reattach the opposite side of shield cap to the tip of a needle tip.
 - After applying, light cure each layer entirely in fanning motion for 30 seconds using a curing unit (halogen lamp) or for 10 seconds using a high intensity visible light with minimum intensity of 1000 mW/cm2. Light curing time may differ depending on shade and depth of the cavity (refer to the above chart for the depth of cure and shade). When the cavity surface is large, light cure in segment.
 - Use dental filling restorative material such as BEAUTIFIL to fill in, light cure, and polish and finish with "Dura White Stones", "T&F Hybrid Points", or "Super-Snap Strips" for interproximal removal and after adjusting occlusion, use polishers such as "CompoMaster", "Super-Snaps" or "OneGloss" for final polishing.

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- For repair of fractured prosthetics/additional filling/CR laminate veneer**
- Prepare according to types of material (resin, porcelain, or tooth).
 - Follow a manufacturer's instructions for the application and curing procedure of etching, priming, and bonding. Fluoride release FL-BOND is recommended.
 - To promote adhesion the use of SHO-FU's FL-Bond is preferably recommended for double action of fluoride release.
 - In any case follow manufacturer's instructions.
 - Pull the shield cap off from a syringe and attach a needle tip by turning it clockwise until it is securely attached. Dispense necessary amount of BEAUTIFIL Flow directly into the cavity. BEAUTIFIL Flow may also be extruded onto a dispensing pad to be applied into cavity with a suitable instrument. Incremental filling of 2mm is recommended. Remove the needle tip from a syringe by turning it counter-clockwise and discard. Then, attach a new needle tip and reattach the opposite side of shield cap to the tip of a needle tip.
 - After applying, light cure each layer entirely in fanning motion for 30 seconds using a curing unit (halogen lamp) or for 10 seconds using a high intensity visible light with minimum intensity of 1000 mW/cm2. Light curing time may differ depending on shade and depth of the cavity (refer to the above chart for the depth of cure and shade). When the cavity surface is large, light cure in segment.
 - Remove excess material with suitable finishers such as "Dura White Stones" or "T&F Hybrid Points", and after adjusting occlusion, use polishers such as "CompoMaster", "Super-Snaps" or "OneGloss" for final polishing.

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- Safety recommendations**
- Avoid using this material to a patient known to have allergies to metaacrylate monomer. (This pertains also to the operator.)
 - If allergic reactions occurs such as eruption or skin inflammation on a patient while using this material, immediately discontinue the use and have the patient seek medical advice. (The same precaution is applied for the operators.)
 - Avoid contacting soft tissues, skin or eyes. In case of contact, rinse immediately with copious amounts of water and seek medical advice, if necessary.
 - Wear UV protective glasses when light curing.
 - Do not use this product in any other purpose except the purpose mentioned in the indication. This product is intended for dental use only.

Compositions	
Base resin :	Bis-GMA/TEGDMA resin
Filler :	multi-functional glass fillers, Improved S-PRG fillers based on fluoro-boroaluminosilicate glass
Filler loading :	54.5 wt%, 34.6 vol% (F02) and 53.8 wt%, 33.3 vol% (F10)
Particle range :	0.01 to 2,0 µm
Mean particle size :	0,8 µm

Storage
Store in a dark room at room temperature (1-30 °C/34-86 °F) avoiding from high temperature, high humidity, direct sunlight and any open flames.

Note: Use this product within the expiration date.
Packaging
Syringe: 2,0g

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